

## ***Deep-amplicon sequencing protocol for fecal nematode identification***

*SMack Lab 2019 (adapted from Avramenko et al. 2015)*

Following extraction with the MoBio DNEasy Blood and Tissue Kit (larval DNA) or Powerlyser Powersoil kit (fecal DNA) as per the manufacturer protocols.

### **Reagents:**

- AMPure XP Beads or suitable alternative
- Forward and reverse adapted primers (Pafco et al., 2018): Strongyl\_ITS-2\_F (acg tct ggt tca ggg ttg) and Strongyl\_ITS-2\_R (atg ctt aag ttc agc ggg ta) (diluted to 10 $\mu$ M concentration).
- NEBNext<sup>®</sup> Ultra<sup>™</sup> II Q5<sup>®</sup> Master Mix
- Illumina barcoding primers
  - Illumina P5 sequence with NEBNext adaptor site and 10 bp i5 index
  - Illumina P7 sequence with NEBNext adaptor site and 10 bp i7 index

**PCR Round One.** This phase puts primers on for our target region (ITS-2).

- Make a master mix as follows (19 $\mu$ L per reaction)
  - 10 $\mu$ L NEBnext
  - 1 $\mu$ L Strongyl\_ITS-2\_F adaptor pool (10 $\mu$ M)
  - 1 $\mu$ L Strongyl\_ITS-2\_R adaptor pool (10 $\mu$ M)
  - 7 $\mu$ L nuclease-free H<sub>2</sub>O.
- Add 19 $\mu$ L of the master mix into each well.
- Add 1 $\mu$ L of sample DNA per well.
- Mix and spin
- PCR I protocol: 95°C for 3 minutes (1x), 98°C for 20 seconds, 62°C for 15 seconds, and 72°C for 15 seconds (25x), and then 72°C for 2 minutes.
- Perform a 1x bead-clean with room temperature ampure beads.
  - Add 20 $\mu$ L of beads to each sample, mix thoroughly
  - Incubate at room temperature for 5 minutes
  - Place on magnetic separator, remove and discard liquid
  - Add 200 $\mu$ L of 70% ethanol, mix to wash beads
  - Place on magnetic separator, remove and discard liquid
  - Repeat ethanol wash as above
  - Dry beads at room temperature on magnet

- Elute in 12 $\mu$ L of water, mix thoroughly
- Incubate at room temperature for 2 minutes
- Place on magnetic separator, transfer 10 $\mu$ L of liquid to a new plate.

**PCR Round Two:** This phase adds Illumina adaptors and barcodes.

- Make a master mix as follows (15.5 $\mu$ L per reaction)
  - 10 $\mu$ L Nebnext buffer
  - 5.5 $\mu$ L H<sub>2</sub>O
- Add 15.5 $\mu$ L of master mix into each well.
- Add unique combination of Illumina P7 (1 $\mu$ L) + P5 (1 $\mu$ L) primers into each well (keep close track of unique combinations).
- Add 2 $\mu$ L of sample DNA per well.
- Mix and spin.
- PCR II protocol: 98°C for 45 seconds (1x), followed by 98°C for 20s, 63°C for 20s (7x), 72°C for two minutes (1x).
- Perform a 1x bead-clean with room temperature ampure beads.
  - Add 20 $\mu$ L of beads to each sample, mix thoroughly
  - Incubate at room temperature for 5 minutes
  - Place on magnetic separator, remove and discard liquid
  - Add 200 $\mu$ L of 70% ethanol, mix to wash beads
  - Place on magnetic separator, remove and discard liquid
  - Repeat ethanol wash as above
  - Dry beads at room temperature on magnet
  - Elute in 22 $\mu$ L of water, mix thoroughly
  - Incubate at room temperature for 2 minutes
  - Place on magnetic separator, transfer 20 $\mu$ L of liquid to a new plate.