

Reduced Representation Bisulfite Sequencing Libraries with EZ-96 DNA Methylation-Lightning™ MagPrep

SMack Lab (adapted from Amanda Lea)

Version 1.2, 7/29/2019

Reagents:

- AMPure XP Beads
- NEB Methylated Illumina Indexes Adaptors

***USE ONLY DILUTED adaptors in subsequent protocol steps. Aliquot a small amount of adaptors into a FRESH tube and dilute 1:20 such that the final concentration of diluted adaptor is ~1 μM.**

**** (MAKE SURE this is the NEB Methylated Illumina adaptor, product number E7535S/L)****

- MspI enzyme (NEB: #R0106)
- 1X NEB CutSmart buffer
- Unmethylated phage DNA (Sigma: #D3654-5UN)
- EB buffer
- NEBNext end prep enzyme + buffer
- NEBNext Blunt/TA Ligase Master Mix _ enhancer
- NEB User Enzyme
- EpiMark Hot Start Taq and buffer
- ZymoEZ-96 DNA Methylation-Lightning™ MagPrep

A few notes:

*Do **not** vortex enzymes or mixtures that contain enzymes

*Be very careful with bisulfite converted DNA. **Do not vortex it or freeze thaw it.** The DNA is single stranded, and therefore very fragile.

*Deep well plates will be needed for lightning conversion steps (included in kit)

Step 1: Fragment DNA

1. Prepare the following mix:

Mix	Storage Location	Volume per Sample	Volume for Mix
10x NEB Cut Smart Buffer		2.5 μL	
MspI Enzyme		1 μL	
0.2 ng/ul unmethylated phage DNA		0.5 μL	

1.1. Add 4 μL of mix to the template DNA (template DNA = 300ng template DNA + nuclease free H₂O for a total of 26 μL)

1.2. Incubate the DNA at 37°C overnight with (do not heat the lid higher than 37°C if using the thermocycler)

***Safe Stopping Point: Freeze the digested DNA at -20°C**

Step 2

2. Prepare the following mix:

Mix	Storage Location	Volume per Sample	Volume for Mix
End Prep Enzyme Mix		1.5 μL	
End Repair Reaction Buffer (10x)		3.25 μL	

- 2.1. Add 4.75 μL of mix to each well of the fragmented DNA (DNA volume is $\sim 30 \mu\text{L}$)
- 2.2. Incubate for 30 minutes @ 20°C , 30 minutes @ 65°C and hold @ 4°C (Run “Neb endprep” program)
- 2.3. Add the following mix to the reaction:

****Make sure that adaptors are diluted 1:20 prior to use in a fresh tube!****

Mix	Storage Location	Volume per Sample	Volume for Mix
NEBNext Methylated Adapter for Illumina (1:20 diluted adapter) **		1.25 μL	
Blunt/TA Ligase Master Mix (Ultra II)		7.5 μL	
Ligation Enhancer		0.5 μL	

- 2.4. Incubate for 20 minutes @ 20°C (Run “NEB Adapter Ligase” program which will hold/pause for User Enzyme Incubation)
- 2.5. Take AMPure beads out of fridge to come to room temperature
- 2.6. Add 1.5 μL USER Enzyme to each sample, pipette to mix, and incubate for 15 minutes @ 37°C (Press enter on “NEB Adapter Ligase program” to end pause and continue the program)

Step 3: Bead-Based Cleanup (1.5x)

- 3.1. Add 66 μL AMPure beads to each sample and gently mix by pipetting up and down
- 3.2. Incubate at RT for 5min
- 3.3. Place plate on magnetic stand until the solution is clear (5 mins)
- 3.4. Remove supernatant (**KEEP BEADS**)
- 3.5. While on the magnetic stand, add 200 μL 80% ethanol, incubate for 30 seconds remove ethanol
- 3.6. Repeat wash: while on the magnetic stand, add 200 μL 80% ethanol, incubate for 30 seconds remove ethanol
- 3.7. After removing the ethanol from the second wash, let the beads dry for ~ 5 -10min or until beads are no longer shiny (**do not over dry** the beads, this will reduce yield)
- 3.8. Remove from magnetic stand and add 21 μL of EB buffer, mix well and incubate 5 minutes @ RT
- 3.9. Place tubes back on magnetic stand and incubate for 5 min @ RT
- 3.10. Transfer 20 μL of supernatant to new PCR strip tube or plate

****Qubit a few samples to check that some still have DNA before proceeding with bisulfite conversion**

***SAFE STOPPING POINT: Freeze the adapter-ligated DNA at -20°C ***

Step 4: Bisulfite Conversion (EZ-96 DNA Methylation-Lightning™ MagPrep)

- 4.1. Add 130 μL of Lightning Conversion Reagent to 20 μL of a DNA sample in a Conversion Plate. Pipette mix 5 times.
- 4.2. Seal the plate with film. Transfer the Conversion Plate to a thermal cycler and perform the following steps:
 - i. 98°C for 8 minutes
 1. 54°C for 60 minutes
 2. 4°C storage (for up to 20 hours)
- 4.3. Pre-heat a plate heating element to 55°C .
- 4.4. Add 600 μL of M-Binding Buffer and 10 μL of MagBinding Beads to each well of a Collection Plate.
- 4.5. Transfer the samples from the Conversion Plate into the Collection Plate (containing the M-Binding Buffer and MagBinding Beads). Pipette mix 5 times.

- 4.6. Let plate stand at room temperature for 5 minutes, then transfer plate to a magnetic stand for an additional 5 minutes or until beads pellet and supernatant is cleared.
- 4.7. With the plate on the magnetic stand, remove the supernatant and discard. **Note:** Some beads will adhere to the sides of the well. Remove supernatant slowly to allow these beads to be pulled to the magnet as the liquid level is lowered.
- 4.8. Remove the Collection Plate from the magnetic stand.
- 4.9. Add 200 μL of M-Wash Buffer to the beads. Pipette mix 5 times.
- 4.10. Place the plate on the magnetic stand for 3 minutes or until beads pellet.
- 4.11. Remove and discard supernatant.
- 4.12. Remove the Collection Plate from the magnetic stand.
- 4.13. Add 200 μL of L-Desulphonation Buffer to the beads. Pipette mix 5 times.
- 4.14. Let plate stand at room temperature (20-30°C) for 15 minutes.
- 4.15. After the incubation, replace the plate on the magnetic stand for 3 minutes or until beads pellet. Remove and discard supernatant. **Important: Take time for handling/re-suspension into account for the total incubation time. Adjust time as necessary to ensure that no sample remains in the L-Desulphonation Buffer for more than 20 minutes**
- 4.16. Remove the Collection Plate from the magnetic stand.
- 4.17. Add 200 μL of M-Wash Buffer to the beads. Pipette mix 5 times.
- 4.18. Place the plate on the magnetic stand for 3 minutes or until beads pellet. Discard supernatant.
- 4.19. Repeat M-Wash Buffer wash once.

Important: Remove as much buffer as possible after final wash to aid in the drying of the beads.

- 4.20. Transfer the plate to a heating element at 55°C for 20-30 minutes to dry the beads and remove residual M-Wash Buffer. Note: Beads will change in appearance from glossy black when still wet to a dull brown when fully dry.
- 4.21. Add 21 μL of M-Elution Buffer directly to the dried beads and pipette mix 5-10 times to re-suspend.
- 4.22. Heat the elution at 55°C for 4 minutes then transfer the plate to the magnetic stand for 1 minute or until beads pellet.
- 4.23. Remove the 20 μL of the supernatant and transfer to a clean Elution Plate.

Note: It is OK if **some** beads are removed with the elution.

The DNA is ready for immediate analysis or can be stored at or below -20°C for later use. For long term storage, store at or below -70°C.

***Quantification should be performed on NanoDrop using RNA settings (Factor 40) or other similar spectroscopy methods using RNA absorbance (If yields are high a second conversion round can be performed)*

SAFE STOPPING POINT: Freeze the digested DNA at -20C

Step 5: PCR Amplification (Indexing) *BE SURE TO USE NEBNext INDEXES*

- 5.1. Make PCR Master Mix as follows:

Mix	Storage Location	Volume per Sample	Volume for Mix
5X EpiMark Hot Start Taq Reaction Buffer		5 μL	
10mM dNTPs		0.5 μL	
EpiMark Hot Start Taq (2 units/ μL)		0.125 μL	

- 5.2. Add 5.625 μL of Master Mix to each well of the Adapter Ligated, BS-treated DNA (has approx. 20 μL of DNA).

- 5.3. Add 0.5 μ L of NEBNext Index Primer i7 to each well.
- 5.4. Add 0.5 μ L of NEBNext Index Primer i5 to each well. Ensure all wells are unique combinations.
- 5.5. Run the following PCR Program:
 1. 95°C for 30 seconds
 2. 16 cycles of:
 1. 95°C for 15 seconds
 2. 61°C for 30 seconds
 3. 68°C for 30 seconds
 3. 68°C for 5 mins
 4. Hold at 4°C

Step 6: Final cleanup and quantification

- 6.1. 2:1 bead clean up
 - a. Add 50 μ L AMPure Beads to each 25 μ L reaction (gently pipetting to mix)
- 6.2. Incubate @ RT for 5 mins
- 6.3. Place plate on magnetic stand until the solution is clear (5 mins)
- 6.4. Remove supernatant (**KEEP BEADS**)
- 6.5. While on the magnetic stand, add 200 μ L 80% ethanol, incubate for 30 seconds remove ethanol
- 6.6. Repeat wash: While on the magnetic stand, add 200 μ L 80% ethanol, incubate for 30 seconds remove ethanol
- 6.7. After removing the ethanol from the second wash, let the beads dry for ~5-10 min or until beads are no longer shiny (**do not over dry** the beads, this will reduce yield)
- 6.8. Remove from magnetic stand and add 22 μ L of nuclease-free H₂O, mix very well and incubate 5 minutes @ RT
- 6.9. Place tubes back on magnetic stand and incubate for 5 min @ RT
- 6.10. Transfer 21 μ L of supernatant to new “candy cane” tube or **plate**
- 6.11. Quantify and visualize libraries on Fragment Analyzer